

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **FIX&PERM Kit 1000 Cell Fixation and Permeabilization Kit (CE/IVD) NMB-GAS-002-1**

Artikelname	FIX&PERM Kit 1000 Cell Fixation and Permeabilization Kit (CE/IVD)
Artikelnummer	NMB-GAS-002-1
Hersteller Artikelnummer	GAS-002-1
Alternativnummer	NMB-GAS-002-1
Hersteller	NordicMubio
Kategorie	Sonstiges
Applikation	FC
Produktbeschreibung	Flow cytometric analyses with monoclonal antibodies were so far mainly restricted to cell surface molecules. Intracellular structures such as cytoplasmic or nuclear enzymes, oncoproteins, cytokines, immunoglobulins etc. were largely excluded from suc...
Puffer	This FIX&PERM Cell Fixation and Permeabilization Kit contains 2 reagents: Fixation Medium (Reagent A) and Permeabilization Medium (Reagent B). It is intended for first fixing cells in suspension with Reagent A and then permeabilizing the cell membrane
Anwendungsbeschreibung	Biological fluids (blood, bone marrow, and others) must be collected under sterile conditions. Anticoagulation with EDTA or heparin is recommended. The samples should be stored at room temperature until used. For optimal results, samples should be proces