

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat anti Bovine Serum proteins, conjugated with FITC, Clone: [Polyclonal], Monoclonal NMB-GAB/TSP/FITC

Artikelname	Goat anti Bovine Serum proteins, conjugated with FITC, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-GAB/TSP/FITC
Hersteller Artikelnummer	GAB/TSP/FITC
Alternativnummer	NMB-GAB/TSP/FITC
Hersteller	NordicMubio
Wirt	Goat
Kategorie	Antikörper
Konjugation	FITC
Produktbeschreibung	Tested in immunoelectrophoresis and double radial immunodiffusion against pooled serum precipitation of at least 12 individual proteins have been observed. Direct fluorescent staining of bovine serum proteins at the cellular and subcellular level of ...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
Puffer	FITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added. It is reconstituted by adding 1 ml sterile distilled water, spun down to remove insoluble particles, divided into sm

Quelle	Pooled whole serum and selected protein fractions isolated from pooled bovine serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formel	FITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added.
Anwendungsbeschreibung	ELISA, Immunocytochemistry, Immunohistochemistry (frozen), (In)direct immunofluorescence.