

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat anti Rat IgG1 (subclass specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal NMB-GARA/IGG1/TRITC

| | |
|--------------------------|--|
| Artikelname | Goat anti Rat IgG1 (subclass specific), conjugated with TRITC, Clone: [Polyclonal], Monoclonal |
| Artikelnummer | NMB-GARA/IGG1/TRITC |
| Hersteller Artikelnummer | GARa/IgG1/TRITC |
| Alternativnummer | NMB-GARA/IGG1/TRITC |
| Hersteller | NordicMubio |
| Wirt | Goat |
| Kategorie | Antikörper |
| Konjugation | TRITC |
| Produktbeschreibung | The reactivity of the antiserum is directed to the subclass IgG1. It does not react with other subclasses of IgG, IgG/Fab fragments, IgM and IgA or any non-Ig protein in rat serum, as tested by immunoelectrophoresis and double radial immunodiffusion.... |
| Klonalität | Monoclonal |
| Klon-Bezeichnung | [Polyclonal] |
| Puffer | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, spun do |

| | |
|------------------------|--|
| Quelle | Pools of purified homogenous IgG1 isolated from pooled rat serum. Freund's complete adjuvant is used in the first step of the immunization procedure. |
| Formel | TRITC-coupled purified hyperimmune goat IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Anwendungsbeschreibung | ELISA, Immunocytochemistry, Immunohistochemistry (frozen), (In)direct immunofluorescence. |