

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit anti Bovine Nucleoside phosphorylase, Clone: [Polyclonal], Monoclonal NMB-NE112/PAB

| | |
|--------------------------|--|
| Artikelname | Rabbit anti Bovine Nucleoside phosphorylase, Clone: [Polyclonal], Monoclonal |
| Artikelnummer | NMB-NE112/PAB |
| Hersteller Artikelnummer | NE112/Pab |
| Alternativnummer | NMB-NE112/PAB |
| Hersteller | NordicMubio |
| Wirt | Rabbit |
| Kategorie | Antikörper |
| Produktbeschreibung | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and... |
| Klonalität | Monoclonal |
| Klon-Bezeichnung | [Polyclonal] |
| UniProt | P55859 |
| Puffer | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. no preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitut |

| | |
|------------------------|--|
| Quelle | Nucleoside phosphorylase isolated and purified from calf spleen. Freund's complete adjuvant is used in the first step of the immunization procedure. |
| Formel | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. no preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Anwendungsbeschreibung | Immunoprecipitation. |