

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit anti Crotalus durissus terrificus venom L-Amino acid oxidase, Clone: [Polyclonal], Monoclonal NMB-NE012/PAB**

Artikelname	Rabbit anti Crotalus durissus terrificus venom L-Amino acid oxidase, Clone: [Polyclonal], Monoclonal
Artikelnummer	NMB-NE012/PAB
Hersteller Artikelnummer	NE012/PAb
Alternativnummer	NMB-NE012/PAB
Hersteller	NordicMubio
Wirt	Rabbit
Kategorie	Antikörper
Applikation	ELISA, ICC, IP, WB
Spezies Reaktivität	Snake
Produktbeschreibung	The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...
Klonalität	Monoclonal
Klon-Bezeichnung	[Polyclonal]
Puffer	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted

Quelle	L-Amino acid oxidase isolated and purified from <i>Crotalus durissus terrificus</i> venom. Freund's complete adjuvant is used in the first step of the immunization procedure.
Formel	Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added.
Anwendungsbeschreibung	Immunoprecipitation.