

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit anti Goat IgG (heavy and light chains) (non Human and Mouse), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal NMB-R7AG/IGG(H+L)/PO

| | |
|--------------------------|---|
| Artikelname | Rabbit anti Goat IgG (heavy and light chains) (non Human and Mouse), conjugated with Horseradish peroxidase, Clone: [Polyclonal], HRP, Monoclonal |
| Artikelnummer | NMB-R7AG/IGG(H+L)/PO |
| Hersteller Artikelnummer | R7AG/IgG(H+L)/PO |
| Alternativnummer | NMB-R7AG/IGG(H+L)/PO |
| Hersteller | NordicMubio |
| Wirt | Rabbit |
| Kategorie | Antikörper |
| Konjugation | HRP |
| Produktbeschreibung | The reactivity of the antiserum is directed to the Fc and Fab subunits of the IgG molecule. It includes a certain degree of reactivity with other immunoglobulins via the common Fab portion. It does not react with any non-Ig protein in goat serum, as ... |
| Klonalität | Monoclonal |
| Klon-Bezeichnung | [Polyclonal] |
| Puffer | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. It is reconstituted by adding 1 ml sterile distilled water, |

| | |
|------------------------|---|
| Quelle | Purified normal IgG isolated from pooled goat serum. Freund's complete adjuvant is used in the first step of the immunization procedure. |
| Formel | Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. |
| Anwendungsbeschreibung | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting. |