

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit anti Porcine Amino acid arylamidase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal NMB-NE011/BIO**

|                          |  |
|--------------------------|--|
| Artikelname              | Rabbit anti Porcine Amino acid arylamidase, conjugated with Biotin, Clone: [Polyclonal], Monoclonal  |
| Artikelnummer            | NMB-NE011/BIO  |
| Hersteller Artikelnummer | NE011/Bio  |
| Alternativnummer         | NMB-NE011/BIO  |
| Hersteller               | NordicMubio  |
| Wirt                     | Rabbit   |
| Kategorie                | Antikörper   |
| Konjugation              | Biotin   |
| Produktbeschreibung      | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...  |
| Klonalität               | Monoclonal   |
| Klon-Bezeichnung         | [Polyclonal]   |
| Puffer                   | Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile di |

|                        |  |
|------------------------|--|
| Quelle                 | Amino acid arylamidase isolated and purified from porcine kidney. Freund's complete adjuvant is used in the first step of the immunization procedure.  |
| Formel                 | Biotin-coupled hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Anwendungsbeschreibung | ELISA, Immunocytochemistry, Immunohistochemistry (paraffin), Dot blot, Immunoblotting.   |