

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Nestin Monoclonal Antibody [10C2] (also known as 2C1.3A11) - Neural Stem Cell Marker - Fluorescein (FITC) / 50ul, IgG1, Mouse IMQ-IQ300FITC**

|                          |   |
|--------------------------|---|
| Artikelname              | Nestin Monoclonal Antibody [10C2] (also known as 2C1.3A11) - Neural Stem Cell Marker - Fluorescein (FITC) / 50ul, IgG1, Mouse |
| Artikelnummer            | IMQ-IQ300FITC   |
| Hersteller Artikelnummer | IQ300FITC   |
| Alternativnummer         | IMQ-IQ300FITC-0.05  |
| Hersteller               | ImmuQuest   |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | FC, ICC   |
| Spezies Reaktivität      | Human   |
| Immunogen                | A 150-amino-acid fragment from the cloned human nestin conjugated to GST  |
| Produktbeschreibung      | Mouse Anti-Nestin [10C2] (also known as 2C1.3A11)A : Fluorescein...   |
| Klonalität               | Monoclonal  |
| Konzentration            | 1mg/ml  |
| Klon-Bezeichnung         | [10C2]  |
| Isotyp                   | IgG1  |
| NCBI                     | <a href="#">10763</a>   |

|                        |  |
|------------------------|--|
| Pubmed                 | <a href=https://pubmed.ncbi.nlm.nih.gov/19658181 target=_blank rel=nofollow>19658181   |
| UniProt                | <a href="#">P48681</a>   |
| Puffer                 | Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-LinkA, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide |
| Reinheit               | Protein G purified   |
| Formulierung           | Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide  |
| Target-Kategorie       | Nestin   |
| Antibody Type          | Monoclonal Antibody  |
| Application Verdünnung | Optimal antibody dilution should be determined by titration.   |