

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Cytokeratin 5/8 Antibody, IgG1, Clone: [SPM268], Mouse, Monoclonal NBT-MSM13X-948-P1ABX

|                          |   |
|--------------------------|---|
| Artikelname              | Cytokeratin 5/8 Antibody, IgG1, Clone: [SPM268], Mouse, Monoclonal  |
| Artikelnummer            | NBT-MSM13X-948-P1ABX  |
| Hersteller Artikelnummer | MSM13X-948-P1ABX  |
| Alternativnummer         | NBT-MSM13X-948-P1ABX-100  |
| Hersteller               | NeoBiotechnologies  |
| Wirt                     | Mouse   |
| Kategorie                | Antikörper  |
| Applikation              | FC, IF, IHC, WB   |
| Spezies Reaktivität      | Bovine, Canine, Ferret, Hamster, Human, Mouse, Porcine, Rat, Sheep  |
| Immunogen                | Cytoskeletal preparation from HeLa cells  |
| Produktbeschreibung      | It reacts with keratin 5 (58kDa) and keratin 8 (52.5kDa). Simple epithelia express cytotkeratin 8 in combination with 18. On the other hand, basal cells of stratified epithelia express cytotkeratin 5 paired with 14. This antibody therefore, reacts with... |
| Klonalität               | Monoclonal  |
| Klon-Bezeichnung         | [SPM268]  |
| Molekulargewicht         | 58kDa (CK5), 52.5kDa (CK8)  |
| Isotyp                   | IgG1  |
| NCBI                     | <a href="#">3852</a>  |

|                        |   |
|------------------------|---|
| UniProt                | <a href="#">P13647</a>  |
| Formulierung           | 0.5 ml at 200 g/ml of antibody Purified from Bioreactor Concentrate by Protein A/G chromatography. Supplied in 10mM PBS, pH 7.4 with 0.05% BSA and 0.05% sodium azide. Contact us if you require this antibody in a different format.                         |
| Antibody Type          | Monoclonal Antibody   |
| Anwendungsbeschreibung | Flow Cytometry (1-2ug/million cells), Immunofluorescence (1-2ug/ml), Western Blot (1-2ug/ml), Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 min at RT)(For staining of formalin-fixed tissues, digest sections with trypsin at 1mg/ml PBS, 15 min at |