

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Bovine Serum Albumin 30% Solution, Unconjugated ROC-BSA-30-1000**

Artikelname	Bovine Serum Albumin 30% Solution, Unconjugated
Artikelnummer	ROC-BSA-30-1000
Hersteller Artikelnummer	BSA-30-1000
Alternativnummer	ROC-BSA-30-1000
Hersteller	Rockland Immunochemicals
Kategorie	Zellen/Zellkultur
Konjugation	Unconjugated
Konzentration	30% by dry weight
Puffer	0.85% (w/v) Sodium Chloride
Formulierung	Liquid (sterile filtered)
Application Verdünnung	ELISA: User Optimized, FLISA: User Optimized, IHC: User Optimized, IF Microscopy: User Optimized, IP: User Optimized, WB: User Optimized
Anwendungsbeschreibung	BOVINE SERUM ALBUMIN 30% Solution is suitable for use in protease sensitive assays such as RIA, EIA and nucleic acid hybridization, use as a stabilizing agent for proteins and enzymes, including dilute solutions of antibody, and use as a blocking agent t