

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 555 ASC-KTL0531

Artikelname	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 555
Artikelnummer	ASC-KTL0531
Hersteller Artikelnummer	KTL0531
Alternativnummer	ASC-KTL0531-1
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Konjugation	AbFluor 555
Produktbeschreibung	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 555 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 555
Anwendungsbeschreibung	LinKine(TM) AbFluor(TM) 555 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 555 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.