

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

LinKine(TM) AbFluor(TM) 594 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 594 ASC-KTL0541

Artikelname	LinKine(TM) AbFluor(TM) 594 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 594
Artikelnummer	ASC-KTL0541
Hersteller Artikelnummer	KTL0541
Alternativnummer	ASC-KTL0541-1
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Konjugation	AbFluor 594
Produktbeschreibung	LinKine(TM) AbFluor(TM) 594 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 594 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 594
Anwendungsbeschreibung	LinKine(TM) AbFluor(TM) 594 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 594 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.