

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 647 ASC-KTL0561**

Artikelname	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 647
Artikelnummer	ASC-KTL0561
Hersteller Artikelnummer	KTL0561
Alternativnummer	ASC-KTL0561-1
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Konjugation	AbFluor 647
Produktbeschreibung	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 647 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 647
Anwendungsbeschreibung	LinKine(TM) AbFluor(TM) 647 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 647 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.