

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 680 ASC-KTL0581

Artikelname	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD), AbFluor 680
Artikelnummer	ASC-KTL0581
Hersteller Artikelnummer	KTL0581
Alternativnummer	ASC-KTL0581-1
Hersteller	Abbkine Scientific
Kategorie	Kits/Assays
Konjugation	AbFluor 680
Produktbeschreibung	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 680 conjugates directly from proteins, peptides, and other ligands that contain a free amino group....
Tag	AbFluor(TM) 680
Anwendungsbeschreibung	LinKine(TM) AbFluor(TM) 680 Labeling Kit (Optimized for samples with molecular weight of 6 KD to 20 KD) is designed for preparing AbFluor(TM) 680 conjugates directly from proteins, peptides, and other ligands that contain a free amino group.