

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

**Anti-Human CD3-FITC/CD19-PE Cocktail, Clone: [OKT-3CB19], FITC/PE,
Monoclonal
ELA-E-AB-FC0004**

Artikelname	Anti-Human CD3-FITC/CD19-PE Cocktail, Clone: [OKT-3CB19], FITC/PE, Monoclonal
Artikelnummer	ELA-E-AB-FC0004
Hersteller Artikelnummer	E-AB-FC0004
Alternativnummer	ELA-E-AB-FC0004-100, ELA-E-AB-FC0004-20, ELA-E-AB-FC0004-100X2
Hersteller	Elabscience
Kategorie	Antikörper
Applikation	FC
Spezies Reaktivität	Human
Konjugation	FITC/PE
Klonalität	Monoclonal
Konzentration	5 uL/Test
Klon-Bezeichnung	[OKT-3CB19]
Puffer	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Anwendungsbeschreibung	For whole blood samples, add 5 uL Anti-Human CD3-FITC/CD19-PE Cocktail to 100 uL anticoagulant-treated blood sample. Mix and incubate the sample at 4C in the dark for 30 min. Remove red blood cells with RBC lysis solution following the manufacturers ins