

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

SARS-CoV-2 Variants Total IgG 6-plex Panel (Flow Cytometry Multiplex Bead Assay, Non-premixed) **ABS-FCM-B08R-96TESTS**

Artikelname	SARS-CoV-2 Variants Total IgG 6-plex Panel (Flow Cytometry Multiplex Bead Assay, Non-premixed)
Artikelnummer	ABS-FCM-B08R-96TESTS
Hersteller Artikelnummer	FCM-B08R-96tests
Alternativnummer	ABS-FCM-B08R-96TESTS
Hersteller	AcroBiosystems
Kategorie	Kits/Assays
Spezies Reaktivität	Virus
Produktbeschreibung	Since December 2019, the emergence of severe acute respiratory syndrome coronavirus 2 (SARSCoV2) and its associated disease, COVID19, has caused a devastating pandemic worldwide.As the virus spreads globally, the continuous emergence of new mutant st...
Quelle	Since December 2019, the emergence of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and its associated disease, COVID-19, has caused a devastating pandemic worldwide.As the virus spreads globally, the continuous emergence of new mutant str
Formulierung	Solid
Target-Kategorie	Spike RBD
Anwendungsbeschreibung	The kit is developed for qualitative detection of Anti-SARS-CoV-2 Variants (Alpha, Beta, Gamma, Delta, Omicron) total IgG antibody in human serum.