

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

AP-TNAP/ALPL Polyclonal Antibody, PerCP Conjugated, Rabbit BSS-BS-1535R-PERCP

Artikelname	AP-TNAP/ALPL Polyclonal Antibody, PerCP Conjugated, Rabbit
Artikelnummer	BSS-BS-1535R-PERCP
Hersteller Artikelnummer	BS-1535R-PERCP
Alternativnummer	BSS-BS-1535R-PERCP-100,BSS-BS-1535R-PERCP-20
Hersteller	Bioss
Wirt	Rabbit
Kategorie	Antikörper
Applikation	IF, WB
Spezies Reaktivität	Human, Mouse
Immunogen	52-150/524
Konjugation	PerCP
Produktbeschreibung	Alkaline phosphatase (ALP) removes phosphate groups from the 5 end of DNA and RNA, and from proteins, at high pH. Most mammals have 4 different isozymes: placental, placental like, intestinal and non tissue specific (found in liver, kidney and bone)....
Klonalität	Polyclonal
Konzentration	1ug/ul
NCBI	11647

Puffer	Aqueous buffered solution containing 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Quelle	KLH conjugated synthetic peptide derived from AP-TNAP/ALPL
Reinheit	Purified by Protein A.
Target-Kategorie	AP-TNAP/ALPL
Application Verdünnung	WB(1:300-5000), IF(IHC-P)(1:50-200), IF(IHC-F)(1:50-200), IF(ICC)(1:50-200)