

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) ASC-BMU109-EN

Artikelname	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red)
Artikelnummer	ASC-BMU109-EN
Hersteller Artikelnummer	BMU109-EN
Alternativnummer	ASC-BMU109-EN-100
Hersteller	Abbkine Scientific
Kategorie	Sonstiges
Applikation	Cell Culture
Produktbeschreibung	SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red) contains 0.25% trypsin and 0.06% EDTA, and the pH value is 7.2-7.4. After filtration and sterilization, the digestive solution can be directly used for the digestion of cultured cells or so...
Tag	Trypsin-EDTA Solution
Formulierung	Liquid solution
Proben	Cells and Tissues
Anwendungsbeschreibung	In the process of using SuperKine(TM) Trypsin-EDTA Solution, 0.25% (With Phenol Red), special attention should be paid to avoid bacterial contamination of the digest. The digestion time of trypsin cell digestive solution should not be too long, otherw