

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

SuperKine(TM) Trypsin-EDTA Solution,0.25% (Without Phenol Red) ASC-BMU110-EN

Artikelname	SuperKine(TM) Trypsin-EDTA Solution,0.25% (Without Phenol Red)
Artikelnummer	ASC-BMU110-EN
Hersteller Artikelnummer	BMU110-EN
Alternativnummer	ASC-BMU110-EN-100
Hersteller	Abbkine Scientific
Kategorie	Sonstiges
Applikation	Cell Culture
Produktbeschreibung	SuperKine(TM) Trypsin-EDTA Solution,0.25% (Without Phenol Red), and the pH value is 7.2-7.4. After filtration and sterilization, the digestive solution can be directly used for the digestion of cultured cells or some tissues....
Tag	Trypsin-EDTA Solution
Formulierung	Liquid solution
Proben	Cells and Tissues
Anwendungsbeschreibung	In the process of using SuperKine(TM) Trypsin-EDTA Solution, 0.25% (Without Phenol Red), special attention should be paid to avoid bacterial contamination of the digest. The digestion time of trypsin cell digestive solution should not be too long, oth