

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

ActiveMax Streptavidin muBeads, premium grade (for cells) ABS-MBS-C009-10MG

| | |
|--------------------------|--|
| Artikelname | ActiveMax Streptavidin muBeads, premium grade (for cells) |
| Artikelnummer | ABS-MBS-C009-10MG |
| Hersteller Artikelnummer | MBS-C009-10mg |
| Alternativnummer | ABS-MBS-C009-10MG |
| Hersteller | AcroBiosystems |
| Kategorie | Zellen/Zellkultur |
| Produktbeschreibung | ActiveMax Streptavidin coupled μ beads, premium grade (for cells) are produced under sterile manufacturing conditions (ISO 5), and no animal or humanderived components are used throughout the production process. It is produced under our rigorous qualit... |
| Quelle | ActiveMax Streptavidin coupled- μ beads, premium grade (for cells) are produced under sterile manufacturing conditions (ISO 5), and no animal- or human-derived components are used throughout the production process. It is produced under our rigorous qualit |
| Formulierung | Powder |
| Target-Kategorie | Streptavidin Beads |
| Anwendungsbeschreibung | ActiveMax Streptavidin coupled- μ beads, premium grade (for cells) can be used for coupling of any biotinylated proteins or antibodies, and are ideal for numerous applications as follows: When coupled with biotinylated CAR target proteins, it can be used t |