

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Chlamydia trachomatis Virulence plasmid protein pGP3-D Protein BYT-ORB1477545

Artikelname	Chlamydia trachomatis Virulence plasmid protein pGP3-D Protein
Artikelnummer	BYT-ORB1477545
Hersteller Artikelnummer	orb1477545
Alternativnummer	BYT-ORB1477545-20,BYT-ORB1477545-100,BYT-ORB1477545-1
Hersteller	Biorbyt
Kategorie	Molekularbiologie
Produktbeschreibung	Recombinant Chlamydia trachomatis Virulence plasmid protein pGP3-D...
Molekulargewicht	40.7 kDa
UniProt	<a href="#">P0CE18</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Chlamydia trachomatis
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	GNSGFYLYNTQNCVFADNIKVGQMTEPLKDQQIILGTTSTPVAAKMTASDGIS LTVSNNPSTNASITIGLDAEKAYQLILEKLGDQILGGIADTIVDSTVQDILDKITT DPSLGLLKAFNNFPITNKIQCNGLFTPRNIETLLGGTEIGKFTVTPKSSGSMFLV SADIIASRMEGGVVLALVREGDSKPYAISYGYSSGVPNLCSLRTRIINTGLTPTT YSLRVGGLESGVWVWVVALSNGNDILGITNTSNVSF
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration