

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rabbit Anti IgG (H&L)-AbFluor 680 AST-SA1025

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname | Rabbit Anti IgG (H&L)-AbFluor 680 |
| Artikelnummer | AST-SA1025 |
| Hersteller Artikelnummer | SA1025 |
| Alternativnummer | AST-SA1025-100, AST-SA1025-1 |
| Hersteller | Assay Biotechnology |
| Wirt | Rabbit |
| Kategorie | Antikörper |
| Applikation | FC, IF, WB |
| Spezies Reaktivität | Mouse |
| Konjugation | AbFluor 680 |
| Produktbeschreibung | Assay Biotechnology secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offe... |
| Quelle | Polyclonal, Rabbit IgG |
| Reinheit | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Formulierung | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |

Application Verdünnung

Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR Odyssey equipment would be recommended.