

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Porphyromonas gingivalis Gingipain R2 (rgpB) BYT-ORB1096114

Artikelname	Recombinant Porphyromonas gingivalis Gingipain R2 (rgpB)
Artikelnummer	BYT-ORB1096114
Hersteller Artikelnummer	orb1096114
Alternativnummer	BYT-ORB1096114-20, BYT-ORB1096114-100, BYT-ORB1096114-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	Recombinant Porphyromonas gingivalis Gingipain R2(rgpB)...
Molekulargewicht	57.1 kDa
UniProt	<a href="#">P95493</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Porphyromonas gingivalis (strain ATCC BAA-308 / W83)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder
Sequenz	YTPVEEKENGRMIVIVPKKYEEDIEDFVDWKNQRGLRTEVKVAEDIASPV TAN AIQQFVKQEYEKEGNLDLYVLLVGDHDKDIPAKITPGIKSDQVYGGQIVGNDHYNE VFIGRFSCESKEDLKTQIDRTIHYERNITTEDKWLGQALCIASAEAGGPSADNGE SDIQHENIIANLLTQYGYTKIICKYDPGVTPKNIIDAFNGGISLANYTGHGSETA WGTSHFGTTHVKQLTNSNQLPFIFDVACVNGDFLYNVP

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration