

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Human Metalloproteinase inhibitor 2 (TIMP2) BYT-ORB1096136

| | |
|--------------------------|--|
| Artikelname | Recombinant Human Metalloproteinase inhibitor 2 (TIMP2) |
| Artikelnummer | BYT-ORB1096136 |
| Hersteller Artikelnummer | orb1096136 |
| Alternativnummer | BYT-ORB1096136-20, BYT-ORB1096136-100, BYT-ORB1096136-1 |
| Hersteller | Biorbyt |
| Kategorie | Proteine/Peptide |
| Produktbeschreibung | Recombinant Human Metalloproteinase inhibitor 2(TIMP2)... |
| Molekulargewicht | 25.1 kDa |
| UniProt | P16035 |
| Puffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Quelle | Homo sapiens (Human) |
| Reinheit | Greater than 90% as determined by SDS-PAGE. |
| Formulierung | Liquid or Lyophilized powder |
| Sequenz | CSCSPVHPQQAF CNADVIRAKAVSEKEVDSGNDIYGNPIKRIQYEIKQIKMFK GPEKDIEFIYTAPSSAVCGVSLDVGKKEYLIAGKAEGDGKMHITLCDFIVPWD TLSTTQKKS LNHR YQMGCECKITRCPMIPCYISSPDECLWMDWVTEKNINGH QAKFFACIKRSDGSCAWYRGAAPPKQEF LDIEDP |

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration