

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Bovine coronavirus Non-structural protein 2a (2a) BYT-ORB1096658

Artikelname	Recombinant Bovine coronavirus Non-structural protein 2a (2a)
Artikelnummer	BYT-ORB1096658
Hersteller Artikelnummer	orb1096658
Alternativnummer	BYT-ORB1096658-20, BYT-ORB1096658-100, BYT-ORB1096658-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	Recombinant Bovine coronavirus Non-structural protein 2a(2a)...
Molekulargewicht	34.4 kDa
UniProt	Q8V438
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Bovine coronavirus (strain 98TXSF-110-LUN) (BCoV-LUN) (BCV)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder
Sequenz	MAVAYADKPNHFINFPLTQFQGFVLNYKGLQFQLLDEGVDCIKQTAPHISLAM LDIQPEDYRSVDVAIQEVIDDMHWGEGFQIKFENPHILGRCIVLDVKGVEELH DDLVNYYIRDKGCVDDQSRKWIGHCTIAQLTDAALSIGKGVDFINSMQFNKITI NPSSPARLEIVKLGAEKKDGFYETIASHWMGIRFEYNPPTDKLAMIMGYCCLEV VRKELEEGDLPENDDAWFKLSYHYENNSWFFRHVYRKSS

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration