

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Bat coronavirus HKU3 Spike glycoprotein (S), partial BYT-ORB1096665

Artikelname	Recombinant Bat coronavirus HKU3 Spike glycoprotein (S), partial
Artikelnummer	BYT-ORB1096665
Hersteller Artikelnummer	orb1096665
Alternativnummer	BYT-ORB1096665-20, BYT-ORB1096665-100, BYT-ORB1096665-500, BYT-ORB1096665-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	Recombinant Bat coronavirus HKU3 Spike glycoprotein(S),partial...
Molekulargewicht	24.3 kDa
UniProt	<a href="#">Q3LZX1</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Bat coronavirus HKU3 (BtCoV) (SARS-like coronavirus HKU3)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder
Sequenz	RVSPTQEVIRFPNITNRCPFDKVFNATRFPNVYAWERTKISDCVADYTVLYNST SFSTFKCYGVSPSKLIDLCFTSVYADTFLIRSSEVRQVAPGETGVIADYNYKLPD DFTGCVIAWNTAKHDTGNYYYYRSHRKTCLKPFDLSSDDGNGVYTLSTYDF NPNVPVAYQATRVVLSFELLNAPATVCGPKLSTELVKNQCVNF

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration