

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Bovine coronavirus Non-structural protein 2a (2a) BYT-ORB1096666

Artikelname	Recombinant Bovine coronavirus Non-structural protein 2a (2a)
Artikelnummer	BYT-ORB1096666
Hersteller Artikelnummer	orb1096666
Alternativnummer	BYT-ORB1096666-20, BYT-ORB1096666-100, BYT-ORB1096666-500, BYT-ORB1096666-1
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Produktbeschreibung	Recombinant Bovine coronavirus Non-structural protein 2a(2a)...
Molekulargewicht	33 kDa
UniProt	<a href="#">P18517</a>
Puffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Quelle	Bovine coronavirus (strain Quebec) (BCoV) (BCV)
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Liquid or Lyophilized powder

Sequenz	MAVAYADKPNHFINFPLTQFEGFVLNYKGLQFQLLDEGVDCIKQTAPHISLAML DIQPEDYRSVDVAIQEVIDDMHWGEGFQIKFDNPHILGRCIVLDVKGVEELHD DLVNYIRDKGCVADQSRKWIGHCTIAQLTDAALSIKENVDFINSMQFNYKITIN PSSPARLEIVKLGAEEKDGFYETIVSHWMGSRFEYNPPTDKLAMIMGYCCSEV VRKELEEGDLPENDDDAWFKLSYHYENNSWFFRHVYRKVL
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration