

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse IgM Isotype Control DyLight 650, Clone: PFR-03 **EXB-L6-803-C100**

Artikelname	Mouse IgM Isotype Control DyLight 650, Clone: PFR-03
Artikelnummer	EXB-L6-803-C100
Hersteller Artikelnummer	L6-803-C100
Alternativnummer	EXB-L6-803-C100
Hersteller	EXBIO
Kategorie	Antikörper
Applikation	IHC, WB, ICC, FC
Konjugation	DyLight™ 650
Konzentration	0.1 mg/ml
Klon-Bezeichnung	PFR-03
Isotyp	Mouse IgM
Puffer	Stabilizing phosphate buffered saline (PBS) solution containing 15 mM sodium azide
Lagerung	2°C to 8°C
Application Verdünnung	Negative control: The reagent is intended as an isotype control to establish the amount of non-specific antibody binding. For your particular experiment, use the same concentration of this control antibody as the recommended working concentration of the

#### Anwendungsbeschreibung

Negative control: The reagent is intended as an isotype control to establish the amount of non-specific antibody binding. For your particular experiment, use the same concentration of this control antibody as the recommended working concentration of the antigen-specific antibody. Also, when working with prediluted antibodies, dilute the isotype control to the same concentration as is the concentration of the antigen-specific antibody in the prediluted antibody solution you are using. If under particular experimental conditions the background signal of the isotype control is too high (usually when working concentrations of used antibodies are above 10 µg/ml of incubation mixture), change the conditions of your experiment to reduce the background.