

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

ELISpot Flex: Canine IFN-gamma (ALP)
MAB-3113-2A

| | |
|--------------------------|--|
| Artikelname | ELISpot Flex: Canine IFN-gamma (ALP) |
| Artikelnummer | MAB-3113-2A |
| Hersteller Artikelnummer | 3113-2A |
| Alternativnummer | MAB-3113-2A |
| Hersteller | MABTECH |
| Kategorie | Kits/Assays |
| Applikation | ELISpot |
| Spezies Reaktivität | Canine |
| Produktbeschreibung | Basic ELISpot kit for enumeration of cells secreting IFN-gamma, The canine IFN-gamma ELISpot kit is ideal for users who want the flexibility of setting up their own assay. The kit enables enumeration of cells secreting canine IFN-gamma, and contains ... |
| Detektionsbereich | Streptavidin-ALP |
| NCBI | 403801 |
| Pubmed | https://pubmed.ncbi.nlm.nih.gov/33466701 target=_blank rel=nofollow>33466701 |
| Target-Kategorie | IFN-gamma |
| Anwendungsbeschreibung | Complementary Products: https://www.mabtech.com/products/elispot-substrate-bcipcnpb-plus-alp-3650-10-0 , https://www.mabtech.com/products/elispot-plates-white-3654-wp . Contents: MT13,MT166, biotin,Streptavidin-ALP |