

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Flex: Porcine IFN-gamma (ALP) MAB-3130-2A

Artikelname	ELISpot Flex: Porcine IFN-gamma (ALP)
Artikelnummer	MAB-3130-2A
Hersteller Artikelnummer	3130-2A
Alternativnummer	MAB-3130-2A
Hersteller	MABTECH
Kategorie	Kits/Assays
Applikation	ELISpot
Spezies Reaktivität	Porcine
Konjugation	ALP
Produktbeschreibung	Basic ELISpot kit for enumeration of cells secreting IFN-gamma, The porcine IFN-gamma ELISpot kit is ideal for users who want the flexibility of setting up their own assay. The kit enables enumeration of cells secreting porcine IFN-gamma, and contain...
Detektionsbereich	Streptavidin-ALP
NCBI	<a href="#">396991</a>
Pubmed	<a href="https://pubmed.ncbi.nlm.nih.gov/28241999">https://pubmed.ncbi.nlm.nih.gov/28241999</a> target=_blank rel=nofollow>28241999
Target-Kategorie	IFN-gamma

Anwendungsbeschreibung	Complementary Products: <a href="https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0">https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0</a> , <a href="https://www.mabtech.com/products/elispot-plates-white-3654-wp">https://www.mabtech.com/products/elispot-plates-white-3654-wp</a> . Contents: pIFN $\gamma$ -I,P2C11, biotin ,Streptavidin-ALP
------------------------	---