

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### ELISpot Flex: Human IgG (ALP) MAB-3850-2A

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | ELISpot Flex: Human IgG (ALP)                                                                                                                                                                                                                                 |
| Artikelnummer            | MAB-3850-2A                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | 3850-2A                                                                                                                                                                                                                                                       |
| Alternativnummer         | MAB-3850-2A                                                                                                                                                                                                                                                   |
| Hersteller               | MABTECH                                                                                                                                                                                                                                                       |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Applikation              | ELISpot                                                                                                                                                                                                                                                       |
| Spezies Reaktivität      | Human, Primate                                                                                                                                                                                                                                                |
| Konjugation              | ALP                                                                                                                                                                                                                                                           |
| Produktbeschreibung      | Basic ELISpot for enumeration of IgG secreting B cells, Human IgG ELISpot Flex enables enumeration of B cells secreting total IgG as well as antigen-specific IgG. It contains a matched pair of monoclonal capture and detection antibodies, Streptavidin... |
| Detektionsbereich        | Streptavidin-ALP                                                                                                                                                                                                                                              |
| NCBI                     | <a href="#">3500</a>                                                                                                                                                                                                                                          |
| Pubmed                   | <a href="https://pubmed.ncbi.nlm.nih.gov/24105029">https://pubmed.ncbi.nlm.nih.gov/24105029</a> target=_blank rel=nofollow>24105029                                                                                                                           |
| Target-Kategorie         | IgG                                                                                                                                                                                                                                                           |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | Complementary Products:<br><a href="https://www.mabtech.com/products/elispot-plates-white-3654-wp">https://www.mabtech.com/products/elispot-plates-white-3654-wp</a> ,<br><a href="https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0">https://www.mabtech.com/products/elispot-substrate-bcipnbt-plus-alp-3650-10-0</a> . Contents: MT91/145,MT78/145, biotin,Streptavidin-ALP,R848,Recombinant human IL-2 |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|