

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Celltase(TM) Cell Detachment Solution ABI-17-1001

| | |
|--------------------------|---|
| Artikelname | Celltase(TM) Cell Detachment Solution |
| Artikelnummer | ABI-17-1001 |
| Hersteller Artikelnummer | 17-1001 |
| Alternativnummer | ABI-17-1001-100ML,ABI-17-1001-500ML |
| Hersteller | Abeomics |
| Kategorie | Zellen/Zellkultur |
| Produktbeschreibung | Celltase is a ready-to-use cell detachment solution of made up of proteolytic and collagenolytic enzymes. Celltase is a direct replacement for Trypsin solution and is very useful for the routine cell disassociation from standard tissue culture plasti... |
| Formel | 1X Celltase enzyme solution containing aqueous buffer, EDTA and Phenol Red. Sterile-filtered, suitable for cell culture. |
| Application Verdünnung | Recommended procedure for general dissociation: Proper aseptic techniques should be followed while handling cell lines and reagents including Celltase. 1. Thaw Celltase at room temperature or at 4C overnight. Gently swirl the Celltase bottle for proper |