

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat IgG anti-Horse IgG (H+L)-unconj., MinX none, Polyclonal DNA-SEC-182548

|                          |   |
|--------------------------|---|
| Artikelname              | Goat IgG anti-Horse IgG (H+L)-unconj., MinX none, Polyclonal  |
| Artikelnummer            | DNA-SEC-182548  |
| Hersteller Artikelnummer | SEC-182548  |
| Alternativnummer         | DNA-SEC-182548  |
| Hersteller               | dianova   |
| Wirt                     | Goat  |
| Kategorie                | Antikörper  |
| Applikation              | WB, IHC, ELISA  |
| Spezies Reaktivität      | Equine  |
| Immunogen                | Horse IgG whole molecule  |
| Konjugation              | Unconjugated  |
| Produktbeschreibung      | Anti-Horse IgG Antibody generated in goat detects horse IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well as fungi ... |
| Klonalität               | Polyclonal  |
| Konzentration            | 10.0 mg/mL  |
| Isotyp                   | Ig  |
| Puffer                   | 0.01 M Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2   |

|                        |  |
|------------------------|--|
| Reinheit               | This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by |
| Formel                 | 10 mM NaPO <sub>4</sub> , 150 mM NaCl, pH 7,2, lyophilisate, Azide/BSA free  |
| Target-Kategorie       | Horse  |
| Antibody Type          | Polyclonal Antibody  |
| Application Verdünnung | WB: 1:2,000 - 1:10,000   |
| Anwendungsbeschreibung | Anti-Horse IgG antibody is suitable for ELISA, western blot, and immunohistochemistry, as well as other assays requiring lot-to-lot consistency.   |