

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit IgG anti-Goat IgG (H+L)-ATTO 425, MinX Ms,Rb, ATTO-425, Polyclonal DNA-SEC-182860**

Artikelname	Rabbit IgG anti-Goat IgG (H+L)-ATTO 425, MinX Ms,Rb, ATTO-425, Polyclonal
Artikelnummer	DNA-SEC-182860
Hersteller Artikelnummer	SEC-182860
Alternativnummer	DNA-SEC-182860
Hersteller	dianova
Wirt	Rabbit
Kategorie	Antikörper
Applikation	WB, IF, FLISA
Spezies Reaktivität	Goat
Immunogen	Goat IgG whole molecule
Konjugation	ATTO-425
Produktbeschreibung	Anti-Goat IgG ATTO dye Antibody generated in rabbit detects goat IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well ...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Goat IgG (H&L) Antibody conjugated to ATTO425 was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelec
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Goat
Antibody Type	Polyclonal Antibody
Application Verdünnung	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 2.0, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Anwendungsbeschreibung	Rabbit Anti-Goat IgG (H&L) Antibody conjugated to ATTO425 is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis