

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Rabbit IgG anti-Goat IgG (H+L)-ATTO 532, MinX Ms,Rb, ATTO-532, Polyclonal DNA-SEC-182867**

Artikelname	Rabbit IgG anti-Goat IgG (H+L)-ATTO 532, MinX Ms,Rb, ATTO-532, Polyclonal
Artikelnummer	DNA-SEC-182867
Hersteller Artikelnummer	SEC-182867
Alternativnummer	DNA-SEC-182867
Hersteller	dianova
Wirt	Rabbit
Kategorie	Antikörper
Applikation	WB
Spezies Reaktivität	Goat
Immunogen	Goat IgG whole molecule
Konjugation	ATTO-532
Produktbeschreibung	Anti-Goat IgG ATTO dye Antibody generated in rabbit detects goat IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well ...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Goat IgG (H&L) conjugated to ATTO 532 was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophore
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Goat
Antibody Type	Polyclonal Antibody
Application Verdünnung	WB: >1:10,000
Anwendungsbeschreibung	Anti-Goat IgG (H&L) conjugated to ATTO 532 has been tested by western blot and is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable f