

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Guinea Pig IgG (H+L)-HRPO, MinX Bo,Ck,Go,Hm,Ho,Hu,Ms,Rb,Rt,Sh, Polyclonal DNA-SEC-182899

| | |
|--------------------------|---|
| Artikelname | Goat IgG anti-Guinea Pig IgG (H+L)-HRPO, MinX Bo,Ck,Go,Hm,Ho,Hu,Ms,Rb,Rt,Sh, Polyclonal |
| Artikelnummer | DNA-SEC-182899 |
| Hersteller Artikelnummer | SEC-182899 |
| Alternativnummer | DNA-SEC-182899 |
| Hersteller | dianova |
| Wirt | Goat |
| Kategorie | Antikörper |
| Applikation | ELISA, WB |
| Spezies Reaktivität | Guinea pig |
| Immunogen | Anti-Guinea Pig IgG (H&L) was produced by repeated immunization with Guinea Pig IgG whole molecule in goat. |
| Konjugation | HRP |
| Produktbeschreibung | Anti-Guinea Pig IgG Peroxidase Antibody generated in goat detects guinea pig IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacter... |
| Klonalität | Polyclonal |
| Konzentration | 1.0 mg/mL |

| | |
|------------------------|--|
| Isotyp | Ig |
| Puffer | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 |
| Reinheit | Anti-Guinea Pig IgG (H&L) peroxidase conjugated antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Guinea Pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. As |
| Formel | 20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% Gentamicin |
| Target-Kategorie | Guinea Pig |
| Antibody Type | Polyclonal Antibody |
| Application Verdünnung | WB: 1:1,000 - 1:5,000 |
| Anwendungsbeschreibung | Anti-Guinea Pig IgG HRP has been tested by ELISA and western blot and is suitable for immunoblotting, western blot, dot blot, ELISA, immunoperoxidase electron microscopy and immunohistochemistry as well as other peroxidase-antibody based enzymatic assays |