

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat IgG anti-Human IgG (H+L)-FITC, MinX Bo,Ck,Go,Gp,Hm,Ho,Ms,Rb,Rt,Sh, Polyclonal DNA-SEC-182991

Artikelname	Goat IgG anti-Human IgG (H+L)-FITC, MinX Bo,Ck,Go,Gp,Hm,Ho,Ms,Rb,Rt,Sh, Polyclonal
Artikelnummer	DNA-SEC-182991
Hersteller Artikelnummer	SEC-182991
Alternativnummer	DNA-SEC-182991
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	DOT
Spezies Reaktivität	Human
Immunogen	Human IgG whole molecule
Konjugation	FITC
Produktbeschreibung	Anti-Human IgG (H&L) Fluorescein generated in goat detects human Immunoglobulin G (IgG), both heavy and light chains of the antibody molecule are present. It is a protein complex composed of four peptide chains - two identical heavy chains and two id...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Secondary Antibodies are prepared from monospecific antiserum by immunoaffinity chromatography using Human IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Human
Antibody Type	Polyclonal Antibody
Application Verdünnung	FLISA Dilution: 1:10,000 - 1:50,000, Flow Cytometry Dilution: 1:500 - 1:2,500, Fluorochrome Protein Value: 2.19, IF Microscopy Dilution: 1:1,000 - 1:5,000
Anwendungsbeschreibung	Anti-Human IgG (H&L) Fluorescein Conjugated Antibody has been tested by dot blot. Secondary antibodies aid in the detection or sorting of target antigens by binding to the primary antibody. Anti-Human IgG (H&L) Fluorescein antibody is ideal for Flow Cyto