

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Rabbit IgG (H+L)-Alk. Phos., MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ALP, Polyclonal DNA-SEC-183347

Artikelname	Goat IgG anti-Rabbit IgG (H+L)-Alk. Phos., MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ALP, Polyclonal
Artikelnummer	DNA-SEC-183347
Hersteller Artikelnummer	SEC-183347
Alternativnummer	DNA-SEC-183347
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	DOT, ELISA
Spezies Reaktivität	Rabbit
Immunogen	Rabbit IgG whole molecule
Konjugation	ALP
Produktbeschreibung	Anti-Rabbit IgG (H&L) Alkaline Phosphatase Antibody generated in goat detects reactivity to Rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G bi...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.05 M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol, pH 8.0
Reinheit	Rabbit IgG (H&L) Alkaline Phosphatase Pre-Adsorbed Conjugated secondary antibody was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwan
Formel	50 mM TrisHCl,150 mM NaCl,1 mM MgCl,0,1 mM ZnCl,50% (v/v) Glycerol,pH 8,0,sterile filtered,0,1% NaN3
Target-Kategorie	Rabbit
Antibody Type	Polyclonal Antibody
Application Verdünnung	WB: 1:500 - 1:2,500
Anwendungsbeschreibung	Anti-Rabbit IgG (H&L) Alkaline Phosphatase Antibody has been tested by ELISA and dot blot and is ideal for Western Blotting, Immunohistochemistry, ELISA as well as other antibody detection methods.