

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Goat IgG anti-Rabbit IgG (H+L)-ATTO 550, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-550, Polyclonal DNA-SEC-183385

Artikelname	Goat IgG anti-Rabbit IgG (H+L)-ATTO 550, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-550, Polyclonal
Artikelnummer	DNA-SEC-183385
Hersteller Artikelnummer	SEC-183385
Alternativnummer	DNA-SEC-183385
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	DOT, WB
Spezies Reaktivität	Rabbit
Immunogen	Rabbit IgG whole molecule
Konjugation	ATTO-550
Produktbeschreibung	Anti-Rabbit IgG (H&L) ATTO 550 Antibody generated in goat detects reactivity to Rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to virus...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Rabbit IgG (H&L) Antibody ATTO 550 was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis
Formel	20 mM K <sub>3</sub> PO <sub>4</sub> , 150 mM NaCl, pH 7.2, lyophilisate, 0.01% NaN <sub>3</sub>
Target-Kategorie	Rabbit
Antibody Type	Polyclonal Antibody
Application Verdünnung	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 2.5, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Anwendungsbeschreibung	Anti-Rabbit IgG (H&L) conjugated to ATTO 550 has been tested by dot blot and western blot and is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is