

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Rabbit IgG (H+L)-ATTO 655, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-655, Polyclonal DNA-SEC-183392

Artikelname	Goat IgG anti-Rabbit IgG (H+L)-ATTO 655, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-655, Polyclonal
Artikelnummer	DNA-SEC-183392
Hersteller Artikelnummer	SEC-183392
Alternativnummer	DNA-SEC-183392
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	WB
Spezies Reaktivität	Rabbit
Immunogen	Rabbit IgG whole molecule
Konjugation	ATTO-655
Produktbeschreibung	Anti-Rabbit IgG (H&L) ATTO 655 Antibody generated in goat detects reactivity to Rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to virus...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Rabbit IgG (H&L) Antibody ATTO 655 was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis
Formel	20 mM K ₃ PO ₄ , 150 mM NaCl, pH 7.2, lyophilisate, 0.01% NaN ₃
Target-Kategorie	Rabbit
Antibody Type	Polyclonal Antibody
Application Verdünnung	WB: >1:10,000
Anwendungsbeschreibung	Anti-Rabbit IgG (H&L) conjugated to ATTO 655 has been tested by western blot and is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable