

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Goat IgG anti-Rat IgG (H+L)-ATTO 425, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rb,Sh, ATTO-425, Polyclonal DNA-SEC-183472

Artikelname	Goat IgG anti-Rat IgG (H+L)-ATTO 425, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rb,Sh, ATTO-425, Polyclonal
Artikelnummer	DNA-SEC-183472
Hersteller Artikelnummer	SEC-183472
Alternativnummer	DNA-SEC-183472
Hersteller	dianova
Wirt	Goat
Kategorie	Antikörper
Applikation	DOT
Spezies Reaktivität	Rat
Immunogen	Rat IgG whole molecule
Konjugation	ATTO-425
Produktbeschreibung	Anti-Rat IgG (H&L) conjugated to ATTO 425 is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including ...
Klonalität	Polyclonal
Konzentration	1.0 mg/mL
Isotyp	Ig

Puffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Reinheit	Rat IgG (H&L) Antibody ATTO 425 was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis res
Formel	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target-Kategorie	Rat
Antibody Type	Polyclonal Antibody
Application Verdünnung	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 2.5, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Anwendungsbeschreibung	Anti-Rat IgG (H&L) conjugated to ATTO 425 has been tested by dot blot and is designed for fluorescence and IF microscopy. The emission spectra for this ATTO conjugate matches the principle output wavelengths of most common fluorescence instrumentation.