

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Dangerous goods 8-element ICP-MS Working Level Tuning Solution of ICP-MS-TS-21 for Thermo iCAP Q, Mix 20 at 1 µg/L Contains: HPS-ICP-MS-TS-20-100

| | |
|----------------------------|---|
| Article Name | Dangerous goods 8-element ICP-MS Working Level Tuning Solution of ICP-MS-TS-21 for Thermo iCAP Q, Mix 20 at 1 µg/L Contains: |
| Biozol Catalog Number | HPS-ICP-MS-TS-20-100 |
| Supplier Catalog Number | ICP-MS-TS-20-100 |
| Alternative Catalog Number | HPS-ICP-MS-TS-20-100-100 |
| Manufacturer | High Purity Standards |
| Category | Biochemikalien |
| Product Description | Q Tune Solution in 2% HNO ₃ + 0.5%HCl**, Barium, Bismuth, Cerium, Cobalt, Indium, Lead, Lithium, Uranium in 2% HNO ₃ + 0.5% HCl. 12 months expiry date. Traceable to NIST 31XX series. ISO 9001:2015 certified, ISO/IEC 17025:2017 and ISO 17034:20166 accre... |