



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## CHO-K1 Lysate NBS-PAT-POYPMP-200

Article Name	CHO-K1 Lysate
Biozol Catalog Number	NBS-PAT-POYPMP-200
Supplier Catalog Number	PAT-POYPMP-200
Alternative Catalog Number	NBS-PAT-POYPMP-200
Manufacturer	Nordic BioSite
Category	Zellen/Zellkultur
Application	WB
Species Reactivity	Hamster
Product Description	CHO-K1 cells are derived from the ovary of an adult Chinese hamster, often used in biological and medical research and commercially in the production of therapeutic proteins. They were introduced in the 1960s, are grown as a cultured monolayer and re
Concentration	2 mg/ml
Buffer	Liquid. In SDS-PAGE Sample Buffer (0.12M Tris-HCl pH6.8, 0.4% SDS, 5% Glycerol, 1% beta-mercaptoethanol, 0.02% Bromophenol blue)
Application Notes	Heat lysate to 95C for 5 minutes and rapidly cooling before use.