

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Influenza A H5N1 Hemagglutinin antibody, Clone: [AT2B7], Unconjugated, Mouse, Monoclonal NBS-AAT-6PB3IX-50

| Article Name               | Influenza A H5N1 Hemagglutinin antibody, Clone: [AT2B7],<br>Unconjugated, Mouse, Monoclonal   |
|----------------------------|---|
| Biozol Catalog Number      | NBS-AAT-6PB3IX-50   |
| Supplier Catalog Number    | AAT-6PB3IX-50   |
| Alternative Catalog Number | NBS-AAT-6PB3IX-50   |
| Manufacturer               | Nordic BioSite  |
| Host                       | Mouse   |
| Category                   | Antikörper  |
| Application                | ELISA, WB   |
| Species Reactivity         | Human, Virus  |
| Immunogen                  | Recombinant Influenza A H5N1/HA1 (17-338aa) purified from Baculovirus   |
| Conjugation                | Unconjugated  |
| Product Description        | H5N1 is a subtype of the species Influenza A virus of the Influenzavirus A genus of the Orthomyxoviridae family. consists of single-stranded eight-segment negative-sense genomic RNAs, helical viral ribonucleoprotein (RNP) complexes (RNA segments NP, |
| Clonality                  | Monoclonal  |
| Concentration              | 1 mg/ml   |
| Clone Designation          | [AT2B7]   |

| Isotype           | IgG1 k   |
|-------------------|--|
| NCBI              | 19417  |
| UniProt           | A9YU04   |
| Buffer            | Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.   |
| Source            | mouse  |
| Antibody Type     | Primary Antibody   |
| Application Notes | The antibody has been tested by ELISA and Western blot analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results. |