

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### TE-2 Lysate NBS-PAT-ROPZ6R-200

Article Name	TE-2 Lysate
Biozol Catalog Number	NBS-PAT-ROPZ6R-200
Supplier Catalog Number	PAT-ROPZ6R-200
Alternative Catalog Number	NBS-PAT-ROPZ6R-200
Manufacturer	Nordic BioSite
Category	Zellen/Zellkultur
Application	WB
Species Reactivity	Human
Product Description	TE-2 cells are derived from human squamous cell carcinoma of the esophagus. The TE cell line series are established in late 1970s and early 1980s by Nishihira et al. in Japan, is one of the first esophageal cancer cell line series that was used throu...
Concentration	2 mg/ml
Buffer	Liquid. In SDS-PAGE Sample Buffer (0.12M Tris-HCl pH6.8, 0.4% SDS, 5% Glycerol, 1% beta-mercaptoethanol, 0.02% Bromophenol blue)
Application Notes	Heat lysate to 95C for 5 minutes and rapidly cooling before use.