

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

kidney (MDBK) cells, Unconjugated NMB-X1007

| | |
|----------------------------|---|
| Article Name | kidney (MDBK) cells, Unconjugated |
| Biozol Catalog Number | NMB-X1007 |
| Supplier Catalog Number | X1007 |
| Alternative Catalog Number | NMB-X1007 |
| Manufacturer | NordicMubio |
| Category | Proteine/Peptide |
| Application | WB |
| Conjugation | Unconjugated |
| Product Description | Total cell lysates are useful as both positive and negative controls in immunoblotting. A431, a human epidermoid carcinoma cell line, and the EGF-stimulated A431 lysates are used as negative and positive controls, respectively, when studying the phos... |
| Concentration | See vial for concentration |
| Buffer | Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel. |
| Formula | Provided as solution in Laemmli electrophoresis sample buffer. Samples need to be boiled for 3-5 minutes, before loading the gel. |
| Application Notes | Suggested use: A volume of 10-25ul is suggested for immunoblotting. |