



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Bovine Serum albumin NMB-BALB**

Article Name	Bovine Serum albumin
Biozol Catalog Number	NMB-BALB
Supplier Catalog Number	BAIb
Alternative Catalog Number	NMB-BALB
Manufacturer	NordicMubio
Host	Bovine
Category	Proteine/Peptide
Product Description	Purified bovine albumin may be used as a reference antigen, calibrator, coating protein and blocking agent in a variety of immunoassays including immunodiffusion, immunoelectrophoresis, ELISA, Western blotting, dot-immunobinding assays (DIBA), haemag
Buffer	Lyophilized, stable white powder. After opening the vial, the lyophilized content is reconstituted by adding 2 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained.
Formula	Lyophilized, stable white powder.