

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Bovine Serum albumin Fr.V, conjugated with FITC NMB-BSA/FITC

Article Name	Bovine Serum albumin Fr.V, conjugated with FITC
Biozol Catalog Number	NMB-BSA/FITC
Supplier Catalog Number	BSA/FITC
Alternative Catalog Number	NMB-BSA/FITC
Manufacturer	NordicMubio
Host	Bovine
Category	Proteine/Peptide
Conjugation	FITC
Product Description	As an efficient counterstain for cell an tissue substrates with immunoconjugates with a label contrasting with the greenish-yellow fluorescence of BSA/FITC. The purpose of counterstaining is to enhance the optical contrast of the immunospecific stain
Buffer	Fluorochrome-coupled bovine serum albumin, fraction V, lyophilized from a solution in phosphate buffered saline (PBS pH 7.2). It is reconstituted by adding 1 ml sterile distilled water.
Formula	Fluorochrome-coupled bovine serum albumin, fraction V, lyophilized from a solution in phosphate buffered saline (PBS pH 7.2).
Application Notes	ELISA,Immunocytochemistry,Immunohistochemistry (frozen),(In)direct immunofluorescence.