

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Chicken anti Dog IgG (H+L), conjugated with FITC, Fluorescein, Gallus NMB-ADIG-F

Article Name	Chicken anti Dog IgG (H+L), conjugated with FITC, Fluorescein, Gallus
Biozol Catalog Number	NMB-ADIG-F
Supplier Catalog Number	ADIG-F
Alternative Catalog Number	NMB-ADIG-F
Manufacturer	NordicMubio
Host	Gallus
Category	Antikörper
Application	IHC
Conjugation	Fluorescein
Product Description	Chicken anti Dog IgG (H+L), conjugated with FITC, Gallus 0,5 mg
Concentration	See vial for concentration.
Buffer	Phosphate buffered saline, pH 7.2, 50% glycerol with 0.075% sodium azide
Source	Chickens were immunized with highly purified dog IgG isolated from normal dog serum. After multiple immunizations in Freunds adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immun
Purity	Egg-yolk derived IgY Affinity Purified
Formula	Phosphate buffered saline, pH 7.2, 50% glycerol with 0.075% sodium azide

۸	nn	licatio	n Ni	otoc
$\sim$	νν	licatio	11 17	Oles

Optimal conjugate working dilutions should be determined for your particular assay conditions.