

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Chicken anti Keyhole Limpet Hemocyanin, conjugated with FITC, Fluorescein, Gallus NMB-AKLH-F

Article Name	Chicken anti Keyhole Limpet Hemocyanin, conjugated with FITC, Fluorescein, Gallus
Biozol Catalog Number	NMB-AKLH-F
Supplier Catalog Number	AKLH-F
Alternative Catalog Number	NMB-AKLH-F
Manufacturer	NordicMubio
Host	Gallus
Category	Antikörper
Application	IHC
Conjugation	Fluorescein
Product Description	Chicken anti Keyhole Limpet Hemocyanin, conjugated with FITC, Gallus 0,5 mg
Concentration	See vial for concentration.
Buffer	Phosphate buffered saline, pH 7.2, 50% glycerol with 0.075% sodium azide
Source	Chickens were immunized with purified Keyhole Limpet Hemocynanin (KLH). After multiple immunizations in Freunds adjuvant, eggs were collected and IgY purified. Anti-KLH was immunoaffinity purified using purified KLH immobilized on a solid support. This a
Purity	Egg-yolk derived IgY Affinity Purified

Formula	Phosphate buffered saline, pH 7.2, 50% glycerol with 0.075% sodium azide
Application Notes	Optimal working dilutions should be determined for your particular assay conditions. Molar F/P ratio is between 3 and 6.